

Automated Filter Integrity Tester

For Radiopharmaceutical Quality Control



www.iphase.com.au

Automated Filter Integrity Tester



The Filter Integrity Tester is a compact standalone bench-top unit capable of fully automating the process of bubble-point testing commonly used vented & non-vented sterilizing filters used in radiopharmaceutical quality control.

The Filter Integrity Tester is designed to **replace manual bench top systems or manual methods of testing** sterilizing filters to enable the following:

- ✓ Automate testing with calibrated sensors.
- ✓ Eliminate operator subjectivity of measurement.
- ✓ Minimise operators exposure to radioactivity.
- ✓ Record & store results in non-editable log.

Key Benefits

Calibrated pressure sensors are used to measure gas pressure and deliver the **true bubble-point**.

Filters are wetted prior to the test to standardize the test to the filter manufacturers specifications which eliminates the need for any product specific bubble-point pressure calibration.

Installation independent results

No need to calibrate pressure drop between the test devices pressure meter and the filter if the filter is located some distance away from the test device.

Vented & Non-Vented Filters

The Filter Integrity Tester is compatible with both vented and non-vented filters as the device automatically detects the filter type and applies the applicable gas flow profile.

This means that from the users perspective both types of filters are tested in the exact same way.

No need to seal vent holes!



Vented filters



Non vented filters

3 simple steps

Step 1



Step 2



Step 1: Select parameter

Select pressure hold test time and bubble-point pressure as per filter manufacturer specification's along with Test ID.

Step 2a: Wet filter

Load up to 35mL of wetting solution (usually water) into the solution reservoir. The solution will be pushed through the filter to wash off residual product which could alter the bubble point.

This step standardises testing and performs the process as per the filter manufactures recommended test procedure.

Step 2b: Pressure hold test


The test is started by pressure being slowly increased up to 80% of the bubble point pressure for a specific amount of time in which the system checks that there are no bubbles/pressure increase on the outlet of the filter.

Step 2c: Bubble Point Test

After the pressure hold test the Filter Integrity Tester slowly increases the pressure to check the bubble point. The Filter Integrity Tester will then record the pressure when bubbles appear/pressure increases on the filter outlet to determine the bubble point pressure.

Step 3: Results

Once the test is complete the filter integrity will display a pass or fail screen with the pressure result.



Simple Installation & Data Logging



Simply connect compressed air, power and waste bottle and its ready to go.

The Filter Integrity Tester is a non-shielded device, if testing radioactive filters then it needs to be located behind sufficient shielding.

Compact dimensions

23cm (H) x 17cm (L) x 17cm (W)

Sensor Calibration



The Filter Integrity Tester requires the two internal pressure sensors to be calibrated once per year. This can be done in one of the two following ways.


Option 1:

Attach a calibrated pressure meter to the device and calibrate the 2 internal pressure sensors using the calibration wizard (this takes about 5-10 minutes)

Option 2:

Ship the Filter Integrity Tester back to iPHASE and we will perform the calibration and have it shipped back to you in no time.

Test History & Data logging



Each test is logged with all details including Date, Time, Test ID, Bubble-point test pressure, Measured bubble-point and Test result.

Data is stored on the device and can be downloaded as a file onto a computer by USB connection.

Technical Specifications

Specifications	
Test Method	Non-destructive true bubble-point using calibrated pressure sensors
Sterile Filters	Non-vented and vented up to 35mm diameter
Filter Wetting	Yes, by loading internal 40mL reservoir with desired solution volume
Gas Supply	Any gas 5-8 bar (73-116 psi), 1/8" OD push in connection
Power Supply	110 V-240 V, 50/60 Hz
Dimensions & Weight	w 15cm x h 23.5cm d 15cm (w 5.9" x h 9.3" x d 5.9") 1.9kg (4.2 lbs)

iPHASE technologies Pty Ltd
Melbourne, Australia
Tel/Fax: +61 (0)3 9884 7596
info@iphase.com.au
www.iphas.com.au



While all care has been taken to ensure that the information contained in this publication is correct, we accept no responsibility for any inaccuracy and reserve the right to modify this information. Technical specifications are based on standard operating conditions and may be subject to variations.